


# Stat-padz®

## Training Electrodes

## For Simulated Defibrillation with ZOLL® Defibrillators




### Instructions for Use

1. Connect electrodes (8900-5005) to the ZOLL Training Cable (8000-1610).
2. Fasten the Training Cable to the Defibrillation Cable.
3. Attach the **Stat-padz** Training Electrodes to any CPR training manikin, pressing the electrodes firmly to the surface and moving any air pockets to the outer edges.
4. Connect the Defibrillation Cable to the Defibrillator and the ZOLL Simulator.

### IMPORTANT:

- Not for clinical use.
- For training purposes only.
- Training system will not connect with the active portions of ZOLL Defibrillators

OPERATING TEMPERATURE: 0°C to 50°C (32°F to 122°F)  
SHORT TERM STORAGE TEMPERATURE: -30°C to 65°C (-22°F to 149°F)  
LONG TERM STORAGE TEMPERATURE: 0°C to 35°C (32°F to 95°F)



NOT MADE WITH  
NATURAL RUBBER LATEX.

R1645-04 Rev. K

# ZOLL®



ZOLL Medical Corporation  
269 Mill Road  
Chelmsford, MA 01824-4105 USA  
800-348-9011 Made in USA



ZOLL International Holding B.V.  
Newsweg 18 6662 PV ELST  
The Netherlands Tel. +31 (0) 481 366410