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
Venoarterial cardiopulmonary bypass for emergency circulatory support²⁾

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Systematic review of percutaneous cardiopulmonary bypass for cardiac arrest or cardiogenic shock states³⁾

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Author	Population	Method of percutaneous cardiopulmonary bypass	Enrolled (number)	Weaned (%)	Survival to discharge (%)
Aliabadi ⁴⁵	Cardiac arrest	Femoral–femoral route with heparin, oxygenator	2	100.0	100.0
Aiba ⁴⁶	Cardiogenic shock complicating acute myocardial infarction	Femoral–femoral route with heparin, oxygenator	26	34.6	19.2
Baird ⁴⁷	Cardiogenic shock (n= 6) or cardiac arrest (n= 19) complicating acute infarction	Femoral–femoral route with heparin, oxygenator	6, 19	50.0, 21.1	33.3, 15.8
Bartlett ⁴⁸	Cardiogenic shock or cardiac arrest	Jugulo–femoral route with heparin, oxygenator	136	n/a	44.0
Bartlett ⁴⁹	Post-cardiotomy cardiogenic shock	Jugulo–femoral route with heparin, oxygenator	6	16.7	16.7
Behringer ⁵⁰	Cardiac arrest 2' to digoxin overdose	Femoral–femoral route with heparin, oxygenator	1	0	0
Bowen ⁵¹	Post-cardiotomy cardiogenic shock or cardiac arrest	Femoral–femoral route with heparin, oxygenator	23	n/a	43.4
Brown ⁵²	Cardiac arrest	Femoral–femoral route with heparin, oxygenator	2	100.0	100.0
Chen ⁵³	Cardiac arrest post-cardiotomy	Femoral–femoral route with unknown heparin, oxygenator	1	100	100
Chen ⁵⁴	Cardiac arrest	Femoral–femoral route with heparin, oxygenator	57	66.7	31.6
Chen ⁵⁵	Cardiogenic shock due to acute myocarditis	Venoarterial route with heparin-bonded circuits, oxygenator	4	100.0	75.0
Cochran ⁵⁶	Cardiac arrest during coronary catheterization	Femoral–femoral route with heparin, oxygenator	3	66.7	66.7
Cofer ⁵⁷	Cardiogenic shock due to acute myocarditis	Venoarterial route with heparin, oxygenator	3	66.7	66.7
Dalton ⁵⁸	Cardiogenic shock (n= 1) or cardiac arrest (n= 4)	Femoral–femoral route with heparin, oxygenator	1, 4	n/a	0, 75.0
Duncan ⁵⁹	Children with cardiac arrest pre- or post-cardiotomy, or during coronary catheterization	Femoral–femoral route, with heparin, oxygenator	11	90.9	63.6
el-Banayosy ⁶⁰	Cardiogenic shock (n= 2) or cardiac arrest (n= 4)	Femoral–femoral route with heparin, oxygenator	2, 4	100.0, 75.0	50.0, 50.0
Frazier ⁶¹	Cardiogenic shock (n= 5), cardiac arrest (n= 2)	Aortic or femoral route with h intra-arterial transvalvular axial flow pump, heparin, no oxygenator	5, 2	100.0, 50.0	60.0, 0
Frizelle ⁶²	Cardiac arrest with massive pulmonary embolus	Femoral–femoral route with TPA	1	0	0
Fumagalli ⁶³	Cardiogenic shock with acute myocardial infarction	Femoral–femoral route with heparin, oxygenator	1	100.0	100.0

Table 1 (continued)

Author	Population	Method of percutaneous cardiopulmonary bypass	Enrolled (number)	Weaned (%)	Survival to discharge (%)
Gacisch ⁶⁴	Cardiogenic shock with acute myocardial infarction	Femoral–femoral route with heparin, oxygenator	3	0	0
Glassman ⁶⁵	Cardiogenic shock	Femoral–atrial route via trans-septal technique with heparin, no oxygenator	3	33.3	0
Grambow ⁶⁶	Cardiogenic shock (n=23) or cardiac arrest (n=7) in catheterization laboratory	Femoral–femoral route with heparin, oxygenator	23, 7	56.5, 0	26.1, 0
Hartz ⁶⁷	Cardiogenic shock (n=3) or cardiac arrest (n=29)	Femoral–femoral route with heparin, no oxygenator	3, 29	100.0, 24.1	100.0, 3.4
Hata ⁶⁸	Cardiogenic shock	Femoral–femoral, Femoral–atrial or aorto-atrial route with heparin, oxygenator	27	51.9	40.7
Hayashi ⁶	Cardiogenic shock post-cardiotomy	Femoral–femoral route with heparin, oxygenator	9	100.0	66.7
Hill ⁶⁹	Cardiogenic shock (n=38) or cardiac arrest (n=111)	Venoarterial route, with unknown heparin, unknown oxygenator	38, 111	n/a	39.5, 13.5
Holley-Jones ⁷⁰	Out-of-hospital cardiac arrest	Femoral–femoral route with heparin, unknown oxygenator	1	100.0	100.0
Holzer ⁷¹	Out-of-hospital cardiac arrest 2' to verapamil overdose	Femoral–femoral route with heparin, unknown oxygenator	1	100.0	100.0
Iguchi ⁷²	Cardiogenic shock with acute myocardial infarction	Femoral–caval route with heparin, oxygenator	1	100.0	100.0
Ihno ⁷³	Cardiogenic shock with acute myocardial infarction	Femoral–femoral route with heparin, oxygenator	1	100.0	100.0
Jaski ⁷⁴	Cardiogenic shock with acute myocardial infarction	Femoral–femoral or jugulo–femoral route with heparin, oxygenator	10	60.0	40.0
Jurmann ⁷⁵	Cardiogenic shock	Femoral–femoral route with heparin, oxygenator	4	25.0	0
Kato ⁷⁶	Cardiogenic shock (n=6) or cardiac arrest (n=3) with myocarditis	Femoral–femoral route with heparin, oxygenator	6, 3	100.0, 100.0	83.3, 66.7
Kanter ⁷⁷	Post-cardiotomy cardiogenic shock	Femoral–femoral route with heparin, oxygenator	4	100.0	100.0
Kawahito ⁷⁸	Post-cardiotomy with shock (n=9) or cardiac arrest (n=4)	Femoral–femoral route with low-dose heparin, heparin-bonded circuits, oxygenator	9, 4	77.8, 100.0	33.3, 50.0
Kawahito ⁷⁹	Myocarditis with shock (n=5) or cardiac arrest (n=1)	Femoral–femoral route with low-dose heparin or heparin-bonded circuits, oxygenator	5, 1	80.0, 100.0	80.0, 100.0

Table 1 (continued)

Author	Population	Method of percutaneous cardiopulmonary bypass	Enrolled (number)	Weaned (%)	Survival to discharge (%)
Obo ¹⁰⁰	Cardiogenic shock or cardiac arrest with acute myocardial infarction	Femoral–femoral route with heparin-coated circuits, oxygenator	21	85.7	42.9
Ohata ¹⁰¹	Cardiogenic shock	Femoral–femoral with heparin-coated circuits, oxygenator	8	62.5	n/a
Orime ¹⁰²	Cardiogenic shock post-cardiotomy or with acute infarction	Femoral–femoral route with heparin, oxygenator	20	35.0	35.0
Pavlidis ¹⁰³	Cardiac arrest in catheterization lab	Femoral–femoral route, unclear heparinization, unclear oxygenator	1	100.0	100.0
Pennington ¹⁰⁴	Cardiogenic shock	Femoral–femoral route, with heparin, oxygenator	14	50.0	28.6
Phillips ¹⁰⁵	Cardiac arrest	Femoral–femoral route with heparin	12	n/a	41.67
Phillips ¹⁰⁶	Cardiogenic shock or cardiac arrest with acute infarction, in catheterization lab or post-cardiotomy	Femoral–femoral with heparin, oxygenator	73	n/a	45.2
Posner ¹⁰⁷	Pediatric cardiac arrest in emergency room	Venoarterial routes with heparin, oxygenator	2	50.0	50.0
Pragliola ¹⁰⁸	Cardiac arrest in catheterization lab	Femoral–femoral route, unknown heparin, oxygenator	3	66.7	33.3
Raithel ¹⁰⁹	Cardiogenic shock with acute infarction (<i>n</i> = 10) or cardiac arrest in catheterization laboratory (<i>n</i> = 10)	Femoral–femoral route with heparin, oxygenator	10, 10	n/a	0, 30.0
Reedy ¹¹⁰	Cardiogenic shock or cardiac arrest in emergency room, catheterization laboratory or intensive care	Femoral–femoral route with heparin, oxygenator	38	63.2	23.7
Reiss ¹¹¹	Cardiogenic shock with myocarditis	Femoral–femoral with unknown heparinization, oxygenator	5	60.0	40.0
Sasaki ¹¹²	Cardiogenic shock (<i>n</i> = 5) or cardiac arrest (<i>n</i> = 4)	Venoarterial route with heparin, oxygenator	5, 4	40.0, 100.0	20.0, 100.0
Sasako ¹¹³	Cardiogenic shock	Femoral–femoral with heparin, oxygenator	40	15.0	7.5
Satoh ¹¹⁴	Cardiogenic shock with myocardial infarction, cardiomyopathy or post-cardiotomy	Femoral–femoral with heparin, oxygenator	4	50.0	25.0
Schwarz ¹¹⁵	Cardiogenic shock not post-cardiotomy (<i>n</i> = 25) or cardiac arrest (<i>n</i> = 21)	Femoral–femoral route with heparin-coated circuits, oxygenator	25, 21	76.0, 42.9	40.0, 14.3
Shawl ¹¹⁶	Cardiogenic shock (<i>n</i> = 3) or cardiac arrest (<i>n</i> = 4) in catheterization laboratory	Femoral–femoral route with unknown heparin, unknown oxygenator	3, 4	66.7, 75.0	33.3, 75.0

Table 1 (continued)


Author	Population	Method of percutaneous cardiopulmonary bypass	Enrolled (number)	Weaned (%)	Survival to discharge (%)
Shaw ¹¹⁷	Cardiogenic shock (<i>n</i> = 6) or cardiac arrest (<i>n</i> = 2) with acute myocardial infarction	Femoral–femoral route with unknown heparin, oxygenator	6, 2	100.0, 50.0	100.0, 50.0
Sugata ¹¹⁸	Cardiogenic shock with acute myocardial infarction	Femoral–femoral route with unknown heparin, oxygenator	1	100.0	100.0
Taferner ¹¹⁹	Neonate with cardiac arrest	Venoarterial, with unknown heparin, oxygenator	1	0	0
Tanaka ¹²⁰	Cardiogenic shock	Femoral–femoral with unknown heparin, oxygenator	17	35.3	23.5
Wakabayashi ¹²¹	Cardiogenic shock	Femoral–femoral, no heparin, no oxygenator	6	83.3	66.7
Walker ¹²²	Cardiogenic shock	Venoarterial, unknown heparin, unknown oxygenator	2	100.0	100.0
Wampler ¹²³	Cardiogenic shock	Femoral–femoral with heparin, unknown oxygenator	41	36.5	31.7
Wang ¹²⁴	Post-cardiotomy cardiogenic shock (<i>n</i> = 8) or cardiac arrest (<i>n</i> = 10)	Femoral–femoral with heparin-coated stents or heparin, oxygenator	8, 10	37.5, 70.0	25.0, 40.0
Wanner ¹²⁵	Cardiac arrest post-PTCA in catheterization lab	Femoral–femoral route with unknown heparin, unknown oxygenator	1	100.0	100.0
Ward ¹²⁶	Cardiogenic shock with myocarditis in intensive care	Venoarterial with heparin, oxygenator	1	100.0	100.0
Willms ¹²⁷	Cardiogenic shock (<i>n</i> = 13) and cardiac arrest (<i>n</i> = 68)	Venoarterial route with heparin or heparin-bonded circuit, oxygenator	13, 68	n/a	23.1, 25.0
Winton ¹²⁸	Post-cardiopulmonary bypass cardiac arrest	Femoral–femoral route with heparin, unknown oxygenator	15	n/a	46.7
Wittenmyer ¹²⁹	Cardiogenic shock (<i>n</i> = 27) or cardiac arrest (<i>n</i> = 73)	Femoral–femoral route with heparin, Bard system	27, 73	77.8, 49.3	11.1, 23.3
Wu ¹³⁰	Cardiogenic shock with acute myocarditis in emergency room	Femoral–femoral route with heparin-coated circuits, hollow fiber oxygenator	1	100.0	100.0
Yasu ¹³¹	Cardiogenic shock with acute myocarditis	Femoral–femoral route with heparin-bonded circuits oxygenator	2	100.0	100.0

Cardiopulmonary resuscitation with assisted extracorporeal life-support versus conventional cardiopulmonary resuscitation in adults with in-hospital cardiac arrest: an observational study and propensity analysis⁴⁾

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Emergency physician-initiated extracorporeal cardiopulmonary resuscitation⁵⁾

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
Table 3
Characteristics of patients who survived ED ECPR with intact neurologic function.

Patient #	Cause of arrest	ECMO duration	Definitive therapy
1	LAD occlusion + refractory VF	43 h 55 min	PCI of LAD + stent
2	RCA occlusion + refractory VF	4 h 25 min	PCI of RCA + stent
3	Severe cardiomyopathy	2 h 3 min	LVAD
4	Hypothermic arrest + refractory VF	24 h 25 min	Re-warming + medical management
5	Aortic dissection, Type A	1 h 15 min	Ascending aortic graft + AV repair


VF, ventricular fibrillation; LAD, left anterior descending coronary artery; RCA, right coronary artery; PCI, percutaneous coronary intervention; LVAD, left ventricular assist device; AV, aortic valve.


Resuscitation of non-postcardiotomy cardiogenic shock or cardiac arrest with extracorporeal life support: The role of bridging to intervention⁶⁾

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Extracorporeal life support associated with hypothermia and normoxemia in refractory cardiac arrest⁷⁾

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
Safety and Feasibility of Prehospital Extra Corporeal Life Support Implementation by Non-Surgeons for Out-of-Hospital Refractory Cardiac Arrest⁸⁾

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
Extracorporeal life support as rescue strategy for out-of-hospital and emergency department cardiac arrest⁹⁾

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Resuscitative extracorporeal membrane oxygenation for in hospital cardiac arrest: A Canadian observational experience¹⁰⁾

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

Refractory cardiac arrest treated with mechanical CPR, hypothermia, ECMO and early reperfusion (the CHEER trial)¹⁾

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
Percutaneous Extracorporeal Life Support in Patients With Circulatory Failure: Results of the German Lifebridge Registry¹¹⁾

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Table 2. Survival at 30 days and procedural endpoints.			
	All [n = 54]	Emergency [n = 39]	Elective [n = 15]
Outcome			
Successfully weaned from support	48 [89%]	33 [85%]	15 [100%]
Survival at 7 days	37 [69%]	22 [56%]	15 [100%]
Survival at 30 days	34 [63%]	19 [49%]	15 [100%]
Procedural endpoints			
Duration of support	162 min [20 min - 26 h]	308 min [20 min - 26 h]	107 min [20 min - 7 h]
Complications			
Hemolysis	0 [0%]	0 [0%]	0 [0%]
Myocardial infarction	1 [2%]	1 [3%]	0 [0%]
SIRS/multi-organ failure	5 [9%]	5 [13%]	0 [0%]
Neurological deficit	2 [4%]	2 [4%]	0 [0%]
Bleeding	2 [4%]	2 [5%]	0 [0%]
Vascular injury	2 [4%]	1 [3%]	1 [7%]
Data presented as n [%] or mean [range].			

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2 LITERATURVERZEICHNIS

1. Stub C, et al. Refractory cardiac arrest treated with mechanical CPR, hypothermia, ECMO and early reperfusion (the CHEER trial). *Resuscitation* 86 (2015) 88-94
<http://dx.doi.org/10.1016/j.resuscitation.2014.09.010>
2. Schwarz B, et al. Venoarterial cardiopulmonary bypass for emergency circulatory support. *Critical Care Medicine* 2003; 31:758-764
3. Nichol G, et al. Systematic review of percutaneous cardiopulmonary bypass for cardiac arrest or cardiogenic shock states. *Resuscitation* 70 (2006) 381-394
<http://dx.doi.org/10.1016/j.resuscitation.2006.01.018>
4. Chen Yih-Sharn, et al. Cardiopulmonary resuscitation with assisted extracorporeal life-support versus conventional cardiopulmonary resuscitation in adults with in-hospital cardiac arrest: an observational study and propensity analysis. *Lancet* 2008;372: 554-61
[http://dx.doi.org/10.1016/S0140-6736\(08\)60958-7](http://dx.doi.org/10.1016/S0140-6736(08)60958-7)
5. Belezio J, et al. Emergency physician-initiated extracorporeal cardiopulmonary resuscitation. *Resuscitation* 83 (2012) 966-970
<http://dx.doi.org/10.1016/j.resuscitation.2012.01.027>
6. Wu M, et al. Resuscitation of non-postcardiotomy cardiogenic shock or cardiac arrest with extracorporeal life support: The role of bridging to intervention. *Resuscitation* 83 (2012) 976-9981
<http://dx.doi.org/10.1016/j.resuscitation.2012.01.010>
7. Fangnoul D, et al. Extracorporeal life support associated with hypothermia and normo-xemia in refractory cardiac arrest. *Resuscitation* 84 (2013) 1519-24
<http://dx.doi.org/doi:10.1016/j.resuscitation.2013.06.016>
8. Lamhaut L, et al. Safety and Feasibility of Prehospital Extra Corporeal Life Support Implementation by Non-Surgeons for Out-of-Hospital Refractory Cardiac Arrest. *Resuscitation* 84 (2013) 1525-9
<http://dx.doi.org/doi:10.1016/j.resuscitation.2013.06.003>
9. Johnson N, et al. Extracorporeal life support as rescue strategy for out-of-hospital and emergency department cardiac arrest. *Resuscitation* 85 (2014) 1527–1532
<http://dx.doi.org/10.1016/j.resuscitation.2014.08.028>
10. Bednarczyk J, et al. Resuscitative extracorporeal membrane oxygenation for in hospital cardiac arrest: A Canadian observational experience. *Resuscitation* 85 (2014) 1713-19
<http://dx.doi.org/10.1016/j.resuscitation.2014.09.026>
11. Jung C, et al. Percutaneous Extracorporeal Life Support in Patients With Circulatory Failure: Results of the German Lifebridge Registry. *J INVASIVE CARDIOL* 2015;27(2): Article in Press

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Copyright 2015
ZOLL Medical Deutschland GmbH
ZOLL Medical Österreich GmbH
ZOLL Medical Switzerland AG