

# Potential Chest Compression Injury Mechanism from Mechanical CPR Systems: Comparison of Load-Distributing Band versus Piston-Driven Systems

TECHNICAL REPORT

NO. 2

It is widely recognized that patient injury can occur during the delivery of chest compressions.<sup>1</sup> Injuries to the ribs, sternum, liver, and spleen have been reported.<sup>2,3,4,5</sup> A review of the injury-related literature for mechanical CPR technologies (Table 1) reveals a higher rate of injury from piston-driven systems for both the overall rate of injury and the difference from manual compressions.

This difference raises the question as to how the respective mechanisms might influence injury rates. Research specific to CPR-related chest compressions demonstrates that surface pressures less than 5.65 pounds per square inch (psi) result in a minimal risk of injury.<sup>6</sup>

**Purpose:** This test was designed to understand the injury-generating potential of different mechanical CPR systems by comparing the forces applied to a chest during compressions.

**Method:** Both piston-driven (LUCAS™ 2) and load-distributing band (AutoPulse®) mechanical CPR systems were applied to a test mannequin in accordance with the manufacturer’s recommendations.<sup>7,8</sup> The anatomical characteristics (chest circumference and stiffness) of the mannequin were calibrated to represent the 50th percentile patient.

Data related to the force of compression at the point of maximal depth was obtained and plotted (Figure 1) with a force mat (Tekscan Body Pressure Measurement System, version 5.83C - Boston, MA) that was placed between the mannequin and the contact surface of the respective CPR systems. Compression depths were recorded with a string potentiometer and plotted using a Dewetron Model No. 3200 (Vienna, Austria). The Compression Force was calculated as follows: maximum depth of compression (inches) x chest stiffness (pounds/inch). The length and width of the Contact Surface Area for each system was measured to the nearest 0.25 inch. The pressure delivered to the chest was determined by dividing the Compression Force by the Contact Surface Area for each system.




Figure 1 – Force plot data for piston-driven (A) and load-distributing band (B) systems.



Figure 2 – The piston-driven (left) and load-distributing band (right) systems shown applied to the calibrated mannequin.

**Table 1 – Reported injury profiles of mechanical CPR technologies\***

	N	CPR Injury Rate	Difference from Manual Compressions
<b>Load-Distributing Band</b> Paradis <sup>9</sup>	1,020	0.8%	Not Reported
<b>Piston Driven</b>			
Wagner <sup>10†</sup>	18	100%	No Control
Smekal <sup>11</sup>	85	58%	+13%
Menzies <sup>12</sup>	57	63%	+ 31%
Obersladstaetter <sup>13</sup>	16	69%	No Control
Englund <sup>14</sup>	221	65% sternal 93% ribs	+132% sternal + 79% rib

\* Excludes superficial injuries such as sternal lacerations, bruising, and abrasions.


† Cath lab experience

**Findings:** Based on an analysis of the data collected one finds:

- The Compression Force is 2.52 times greater for the load-distributing band system than for the piston-driven system; 280 lbs. versus 111 lbs., respectively. (Table 2)
- The Contact Surface Area is 23 times greater for the load-distributing band system than for the piston-driven system; 100 versus 4.4 inch<sup>2</sup>, respectively (Table 2).
- The pressure delivered to the chest at the point of maximum compression is 2.80 psi and 25.45 psi, respectively, for the load-distributing band and piston-driven systems (Figure 3).
  - The pressure delivered by the load-distributing band system is 50% lower than the reported threshold for chest compression injuries.
  - The piston-driven system operates 4.50 times above the reported threshold for safe chest compressions.

**Table 2 – Measure forces and contact surfaces**

Technology	Compression Force (lbs)	Contact Surface Units (in <sup>2</sup> )
Load-Distributing Band	280.0	100.0
Piston-Driven	111.3	4.4



**Figure 3 – Pressure delivered to the chest**

**Conclusion:** The smaller contact surface of the piston-driven system (LUCAS 2) results in 1) a pressure delivered to the chest that substantially exceeds the threshold for compression injury, and 2) the larger contact surface of the load-distributing system (AutoPulse) is designed to provide a generous safety margin as the pressure delivered is well below the threshold for chest compression injuries.

<sup>1</sup> American Heart Association. 2010 Guidelines for CPR and ECC. *Circulation*. 2010;S639-S946.

<sup>2</sup> Clark DT. Complications following closed-chest cardiac massage. *JAMA*. 1962;18:1337-338.

<sup>3</sup> Patterson RH, Burns WA, Janotta FS. Complications of external cardiac resuscitation: a retrospective review and survey of the literature. *Med Ann DC*. 1974;43:389-394.

<sup>4</sup> Nagel EL, Fine EG, Krischer JP, Davis JM. Complications of CPR. *Critical Care Medicine*. 1981;9:424.

<sup>5</sup> Krischer JP, Fine EG, Davis JM, Nagel EL. Complications of cardiac resuscitation. *Chest*. 1987;92:287-291.

<sup>6</sup> Knoell CK. Thoracic response to blunt frontal loading. SH. *Published by Society of Automotive Engineers, Inc.* 1976.

<sup>7</sup> Jolife. LUCAS™ 2 Chest Compression System Instructions for Use (100666-00E). Lund, Sweden. 2009.

<sup>8</sup> ZOLL Circulation. AutoPulse Resuscitation System Model 100 User Guide (11440-001 Rev. 3). Sunnyvale, CA. 2009.

<sup>9</sup> Paradis NA, Kamlan D, Ghilarducci D, Palazzolo J. The California AutoPulse quality assurance registry. *Circulation*. 2009;120:S1457.

<sup>10</sup> Wagner H, Terkelsen CJ, Friberg H, Harnek J, Kern K, Lassen JF, Olivecrona GK. Cardiac arrest in the catheterisation laboratory: a 5-year experience using mechanical chest compressions to facilitate PCI during prolonged resuscitation efforts. *Resuscitation*. 2010; 81:383-387.

<sup>11</sup> Smekal D, Johansson J, Huzevka T, Rubertsson S. No difference in autopsy detected injuries in cardiac arrest patients treated with manual chest compressions compared with mechanical compressions with the LUCAS™ device – a pilot study. *Resuscitation*. 2009; 80:1104-1107.

<sup>12</sup> Menzies D, Barton D, Darcy C, Nolan N. Does the LUCAS device increase trauma during CPR. *Resuscitation*. 2008;77S: AS-034.

<sup>13</sup> Obersladstaetter D, Baubin M, Freund M, Rabl W. Thorax injuries after CPR. *Resuscitation*. 77S; S51, AS-039.

<sup>14</sup> Englund E, Silverstople J, Halvarsson BL. Injuries after cardiopulmonary resuscitation: a comparison between LUCAS mechanical CPR and standard CPR. *Resuscitation*. 2008; 77S: S13, AS-036.