
Solex® Catheter Kit - Instructions for Use - Model SL-2593AE/8700-0671-40
Kit de cathéter Solex® - Mode d'emploi - Modèle SL-2593AE/8700-0671-40
Solex® -Katheter-Kit – Gebrauchsanweisung – Modell SL-2593AE/8700-0671-40
Kit de catéter Solex® - Instrucciones de uso - Modelo SL-2593AE/8700-0671-40
Kit de Cateter Solex® - Instruções de Uso - Modelo SL-2593AE/8700-0671-40
Kit de Cateter Solex® - Instruções de Utilização - Modelo SL-2593AE/8700-0671-40
Kit per catetere Solex® - Istruzioni per l'uso - Modello SL-2593AE/8700-0671-40
Solex® -katheterkit - Gebruiksaanwijzing - Model SL-2593AE/8700-0671-40
Zestaw cewnikowy Solex® - Instrukcja eksploatacji - Model SL-2593AE/8700-0671-40
Sada katetru Solex® - Pokyny pro použití - Model SL-2593AE/8700-0671-40
Komplet katetera Solex® - Upute za uporabu - model SL-2593AE/8700-0671-40
Solex® katéterkészlet – Használati utasítás az SL-2593AE/8700-0671-40 jelű típushoz
Комплект для катетеризации Solex® - Инструкция по эксплуатации - Модель SL-2593AE/8700-0671-40
Solex®-katetersats - Bruksanvisning - Modell SL-2593AE/8700-0671-40
Solex®-katetripakkaus - Käyttöohjeet - Malli SL-2593AE/8700-0671-40
Solex® kateterkit - Brugsanvisning - Model SL-2593AE/8700-0671-40
Solex® kateterio rinkinys - Naudojimo instrukcija - Modelis SL-2593AE/8700-0671-40
Solex® katetra komplekts - Lietošanas instrukcijas - Modelis SL-2593AE/8700-0671-40
Κιτ καθετήρα Solex® - Οδηγίες χρήσης - Μοντέλο SL-2593AE/8700-0671-40

English (en).....	2
Français (fr).....	8
Deutsch (de).....	15
Español (es).....	23
Português (ptbr).....	30
Português (pt).....	37
Italiano (it).....	44
Nederlands (nl).....	51
Polski (pl).....	59
Česky (cs).....	66
Hrvatski (hr).....	74
Magyar (hu).....	81
Русский (ru).....	88
Svenska (sv).....	96
Suomi (fi).....	103
Dansk (da).....	110
Lietuviškai (lt).....	117
Latviski (lv).....	124
Ελληνικά (el).....	131

ZOLL Circulation, Inc.
San Jose, CA 95131
USA
+1-408-541-2140

English (en)

Solex® Intravascular Heat Exchange Catheter (Solex Catheter) Kit includes:

Quantity	E.U.
1	Solex Intravascular Heat Exchange Catheter 9.3 French x 26cm Triple Infusion Lumen Extension Line Clamps Radiopaque Polyurethane Applause™ Heparin Coating
2	Guidewires .032" (0.81 mm) x 65cm
1	Vessel Dilator 10.5 F x .038" (3.6 mm x 1.0 mm)
1	Detachable Suture Tab & Clip
1	18ga x 2 ½" (1.3 mm x 6.3 cm) Radiopaque OTN Catheter
1	000 Silk Suture
6	4" x 4" (10 cm x 10 cm) Gauze Sponges
1	No. 11 Surgical Blade w/ long handle
1	5cc Syringes & 22ga x 1½" (0.7 mm x 38 mm) Needles
1	Fenestrated Drape
1	18ga x 2¾" (1.3 mm x 70 mm) Needle

Device Description

The Solex Catheter is a sterile, single-use, flexible catheter designed for placement in the superior vena cava from an insertion site in the jugular vein. The Solex Catheter serpentine balloon provides a heat exchange surface for the cooling or rewarming of the patient's blood. The serpentine balloon has a cross sectional area of 54mm², which is approximately equal to the cross sectional area of an 8mm OD standard cylindrical balloon. The Solex Catheter is intended for connection to a single-use, disposable CoolGard 3000®/Thermogard XP® Start-Up Kit (supplied separately) and the CoolGard 3000 or Thermogard XP System. An YSI-400 Temperature probe is required for the operation of the CoolGard 3000 or Thermogard XP. A dilator and guidewire are required for the percutaneous insertion of the Solex Catheter. Three (3) lumens are available for infusion and sampling.

Infusion Port	Flow Rate	Priming Volume
Proximal Port (blue)	1300 ml/hr	0.3 cc
Medial Port (white)	1300 ml/hr	0.3 cc
Guidewire Port (brown)	1900 ml/hr	0.4 cc

The Solex Catheter blood contact surfaces are treated with Applause™ Heparin Coating. Applause is a trademark of SurModics, Inc., registered in the U.S. Patent and trademark office.

Sterility

Ethylene oxide sterilized. The Solex Catheter is supplied sterile for single use only and should not be resterilized. The package should be inspected prior to use to ensure that the sterility barrier has not been compromised.

Storage

Store in a cool, dry place.

Indications for Use

The Solex® Intravascular Heat Exchange Catheter in combination with the CoolGard 3000 or Thermogard XP System permits circulation of temperature controlled sterile saline through a heat exchanger to cool / warm the patient's blood in patients for whom the risks of a central line are warranted. Maximum use period: 48 hours.

Contraindications

1. The risks of the catheter are essentially those of a central venous line. The catheter should not be used in patients for whom central line placement is not indicated.
2. Bleeding diathesis.
3. Active sepsis.
4. Infection or active bleeding at the site of catheter insertion.
5. Patients with no vascular access, or a vascular system which will not accommodate catheter including patients with vena cava filters or other implanted impediments to passage of the catheter.
6. Patients for whom the required temperature monitoring cannot be established.
7. Hypothermia is contraindicated in patients who have hematological diseases that will be made worse with hypothermia e.g. any disease that produces cryoglobulinemia, any hemoglobinopathy in which hemolytic anemia can be precipitated by cold including Sickle Cell Disease or Thalassemia.
8. Not intended for pediatric or neonatal use.

Warnings and Precautions

Central venous catheterization should only be performed by well-trained personnel well-versed in anatomical landmarks and safe technique. Personnel should also have knowledge of potential complications.

USE JUGULAR VEIN APPROACH ONLY

Do not attempt to insert another catheter into the same vessel(s) as the Solex Catheter if its tip will reside in the superior vena cava.

If another catheter is to occupy the same vessel(s) as the Solex Catheter, you MUST

- Place both guidewires at the same time, and
- Insert the Solex Catheter AFTER the other catheter to prevent entrapment of the other catheter within the serpentine balloon of the Solex Catheter.

Entrapment of another catheter within the serpentine balloon of the Solex Catheter should be regarded as a surgical emergency.

SINGLE USE ONLY. Product designed for single use only. Do not resterilize or reuse. Do not reinsert, once removed from patient. Do not alter the catheter in any way. Maximum use period: 2 days.

Potential risks with re-use of a single use device include but are not limited to

- Potentially life threatening infection
- Toxic shock due to degradation of materials
- Increased risk of thrombosis
- Reduced heat exchange power
- Device failures

WARNING: Catheter should be placed via a jugular vein approach only. Do not allow catheter to be placed into right atrium or right ventricle. Placement in right atrium or right ventricle can result in severe patient injury or death.

1. **Caution: The use of an introducer sheath is NOT recommended with the Solex Catheter.**
2. Catheter should be positioned so that the distal tip of catheter is in the superior vena cava above its junction with the right atrium and parallel to the vessel wall. X-ray examination should be used to ensure that the catheter is not in the right atrium or ventricle. The distal tip of the catheter should be positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.
3. Cardiac Tamponade: Placement of indwelling catheters in the right atrium is a practice that may lead to cardiac perforation and tamponade. Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size. The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion. Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures, such as aspiration of air emboli during neurosurgery. Such procedures are nevertheless risk prone and should be closely monitored and controlled.
4. Alcohol and acetone can weaken the structure of the polyurethane material. Care should therefore be taken when infusing drugs containing alcohol or when using alcohol or acetone when performing routine catheter care and maintenance. Alcohol should not be used to declot the catheter.
5. Use of a syringe smaller than 10 ml to irrigate or declot an occluded catheter may cause intraluminal leakage or catheter rupture.
6. If blood is observed within the circuit circulating the sterile saline, stop the procedure.
7. The catheter is coated with Heparin. This may induce or aggravate pre-existing Heparin induced thrombocytopenia (HIT).
8. Possible complications with central venous catheters include: atrial or ventricular perforation, cardiac tamponade, air embolism, catheter embolism, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, hematoma formation, hemorrhage, nerve damage and dysrhythmias.
9. All Luer-Lock connections and covers must be securely tightened to prevent air embolism or fluid or blood loss.
10. Never use excessive force in moving the catheter or guidewire. If resistance is encountered, an x-ray should be performed to identify the reason for the resistance.

11. Passage of the guidewire into the right heart can cause dysrhythmias, right bundle branch block, vessel wall, atrial or ventricular perforation.
12. Use only sterile saline for catheter priming and as the circulating fluid in the catheter.
13. Catheter should be routinely inspected for flow rate, security of dressing, correct catheter position and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed. A potential risk associated with multi-lumen balloon catheters is a lack of therapy including failure to infuse drug therapy through the infusion lumens.
14. Only x-ray examination can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If catheter position has changed, perform an x-ray examination to confirm catheter tip position.
15. For blood sampling, temporarily shut off remaining infusion ports through which solutions are being infused.
16. Use only a 30cc or smaller syringe for blood sampling.
17. Use only ZOLL suture tab and clip provided in kit to prevent catheter damage.
18. Do not infuse into the "IN" (INFLOW) and "OUT" (OUTFLOW) Luer-Lock connections as this will result in lack of therapy.
19. Use care when infusing drugs that may be affected by cool temperatures (as low as 4°C). Mannitol containing solutions are temperature sensitive and must not be delivered through the Solex Catheter except for rapid push of up to a concentration of 20% mannitol solution followed by saline flush. Higher than a 20% concentration of mannitol or drip or infusion pump delivery of mannitol must be done via a separate line.
20. **WARNING: When connecting infusion sets/injection systems to ZOLL Catheters do not exceed 100 psi/689 kPa.**
21. For patients being made hypothermic, the hypothermia itself may exacerbate some disease states. Care should be taken to properly monitor patient homeostasis during hypothermia.
 - Cardiac rhythm disturbances – both bradycardia and ventricular tachyarrhythmia.
 - Clotting and coagulations function. Patients at risk for disturbances of their clotting or coagulation function should be closely monitored during hypothermia.
 - Blood gas and pH analysis. Hypothermia modifies resting pH and PaCO₂. Physicians should be aware that of the effect of temperature upon the result.
 - Prolonged hypothermia depresses the immune response and lung function.

WARNING: INTRALUMINAL LEAKAGE

Intraluminal leakage between the saline lumen and infusion lumens is an uncommon but potential catheter failure. In the event of such a failure, sterile saline from the cooling circuit will be introduced into the patient. Intraluminal leakage will usually be associated with a fluid loss alarm which will stop the System. ALWAYS INVESTIGATE FLUID LEVEL ALARMS. The cooling circuit is a closed loop system – usually fluid loss alarms indicate a breach somewhere in this closed loop. With any fluid loss alarm, check both the integrity of the catheter and the Start-Up Kit (see below).

To check the integrity of the catheter

1. Stop operation of the CoolGard 3000 or Thermogard XP System.
2. Disconnect the Start-Up Kit from the catheter and properly cap both the catheter and Start-Up Kit using an aseptic technique.
3. Fill a sterile 10 ml syringe with sterile saline.
4. Connect the syringe to the INFLOW lumen of the catheter and disconnect the OUTFLOW cap. Infuse the 10 ml of saline – it should flow out the OUTFLOW lumen.
5. Now cap the OUTFLOW lumen and pull 5 cc of vacuum and sustain this for at least 10 seconds. Approximately 4 ml of saline (not blood), should enter the syringe and you should be able to maintain the vacuum.
6. Ease the vacuum and recap the INFLOW lumen.

To check the integrity of the Start-Up Kit

1. Look for obvious leakage.
2. Remove the tubing from the peristaltic (roller) pump and inspect for damage (return it to position if not damaged).
3. Check along the tubing from the pump to the patient for sources of fluid loss.
 - Look for damage to the tubing and/or the presence of air within the tubing.
 - Inspect, and tighten as necessary, each Luer fitting (do not use instruments to tighten Luer fittings).
4. Similarly, check the tubing that returns to the pump from the patient. Examine the saline bag to ensure that it has not been accidentally compromised (for example, the spike may have damaged the bag wall).
5. Trace the tubing from the saline bag back to the pump.

Additional warnings and cautions are located in following instructions.

Materials Required

Quantity	Description
1	Solex Catheter Kit for percutaneous introduction
1	Bag of Sterile Saline (not provided)
1	CoolGard 3000/Thermogard XP Start-Up Kit (provided separately)
1	CoolGard 3000 or Thermogard XP System (provided separately)
1	YSI-400 Temperature Probe (not provided)

Catheter Preparation and Insertion

NOTE: The Solex Catheter has radio-opaque marker bands to assist in identification of the catheter during and after insertion when viewed using X-ray equipment. The tip has one marker band and contains barium sulfate to make it radio-opaque. The other end of the balloon has one marker band.

Catheter Preparation Procedure

1. Using sterile technique remove the Solex Catheter from the kit.

2. **Caution: Do not prime the Solex balloon of the Solex Catheter. The Solex balloon is provided fully deflated. There is minimal air within the Solex balloon. It is not possible to force fluid through the collapsed balloon. DO NOT REMOVE THE PROTECTIVE TUBE THAT COVERS THE BALLOON UNTIL IMMEDIATELY PRIOR TO INSERTION.**
3. Remove caps from the INFLOW and OUTFLOW Luer hubs.
4. **WARNING: Never inject positive pressure into the INFLOW hub with the OUTFLOW Luer cap in place.**
5. Using 5 cc or larger syringe, flush the distal, proximal and medial infusion lumens with sterile saline. Clamp or attach injection caps to the proximal and medial infusion lumens. Leave the distal Luer uncapped for guidewire passage.
6. **WARNING: Do not cut the catheter to alter length.**

Catheter Insertion

1. **WARNING: Use jugular vein approach only.**
2. Place patient in a slight Trendelenburg position as tolerated to reduce the risk of air embolism.
3. Prep and drape puncture site as required.
4. **Caution: Always prime catheter infusion lumens before it is inserted into patient.**
5. Obtain jugular venous access using standard percutaneous techniques. Access should be maintained with a .032" (0.81mm) guidewire. See instructions for guidewires.
6. **Caution: Do not use a guidewire larger than .032" (0.81mm) with the Solex Catheter.**
7. **WARNING: Do not attempt to re-insert a partially or completely withdrawn OTN (over the needle) introducer needle from its introducer catheter.**
8. Holding the guidewire in place, remove introducer catheter. **Caution:** Maintain a firm grip on the guidewire at all times.
9. Enlarge the cutaneous puncture site with cutting edge of scalpel positioned away from the guidewire. **WARNING: Do not cut guidewire.** Use vessel dilator to enlarge site as required. Do not leave vessel dilator in place as an indwelling catheter to minimize risk of possible vessel wall perforation.
10. Remove catheter membrane cover.
11. Thread tip of Solex Catheter over guidewire. Maintain a sufficiently firm grip on the guidewire during catheter insertion. Grasping the catheter near the skin, advance catheter into vein.
12. Using centimeter marks on the catheter as positioning reference points, advance catheter to final indwelling position.
13. Hold catheter at desired depth and remove guidewire. If resistance is encountered when attempting to remove the guidewire after catheter placement, the guidewire may be kinked at the tip of the catheter. If resistance is encountered, withdraw the catheter relative to the guidewire about 2-3 cm and attempt to remove the guidewire. If resistance is again encountered remove the guidewire and catheter simultaneously.

14. **Caution: Do not apply undue force to the guidewire.**
15. Verify that the guidewire is intact upon removal.
16. Check lumen placement by attaching a syringe to the distal infusion Luer hub and aspirate until a free flow of venous blood is observed. Connect infusion Luer to appropriate Luer-Lock line as required. Unused infusion port may be "locked" through injection cap using standard hospital protocol. A slide clamp is provided on the tubing to occlude flow through the infusion lumen during line and injection cap changes. **Caution: To minimize risk of damage to the tubing from excessive pressure, the clamp must be opened prior to infusing through the lumen.**
17. **Caution: Do not clamp or occlude INFLOW or OUTFLOW lines. This can cause line blockage and possible failure.**
18. Secure and dress insertion site and catheter temporarily.
19. Verify catheter tip position by chest x-ray immediately after placement. X-ray exam must show the catheter located in the SVC with the distal end of the catheter parallel to the vena cava wall. If the catheter tip is mal-positioned, reposition and re-verify.
20. Proximal radiopaque marker indicates proximal end of balloons to ensure that balloons reside completely in vessel. If catheter is mal-positioned, reposition and re-verify.
21. Secure catheter to patient. Use juncture hub side wings as primary suture site to minimize the risk of catheter migration.
22. The ZOLL suture tab and clip can also be used as an additional attachment point. Assure that catheter body is secure and does not slide.
23. **Caution: Use only the ZOLL suture tab and clip provided in the kit. Catheter damage may result if other tabs or clips are used.**
24. **Caution: Do not suture directly to the outside diameter of the catheter to minimize the risk of cutting or damaging the catheter or impeding catheter flow.**
25. Dress puncture site per hospital protocol. Maintain the insertion site with regular meticulous redressing using aseptic technique.
26. Record on the patient's chart the indwelling catheter length using the centimeter marks on the catheter shaft as reference. Frequent visual reassessment should be made to ensure that the catheter has not moved.
27. Attach a primed Start-Up Kit to Solex Catheter by connecting the male Luer of the Start-Up Kit to the female INFLOW Luer of the Solex Catheter and the female Luer of the Start-Up Kit to the male OUTFLOW Luer of the Solex Catheter. White "ZOLL" tags are fitted to the INFLOW and OUTFLOW extension tubes to help identify them. Assure that a sufficient amount of sterile saline is present at the ends of the hubs to make an air free connection. **Refer to CoolGard 3000 or Thermogard XP manual for details on System operation.**
28. **WARNING: Failure to connect the Start-Up Kit correctly to the catheter could result in catheter failure. Do not attach Start-Up Kit to the brown, white, or blue Infusion Luers.**
29. **Caution: Do not place any stopcocks in line that may be inadvertently shut off. This can cause line blockage and possible failure.**

30. Pump saline through Start-Up Kit and catheter to assure that all connections are secure and that there is no leaking. Allow any remaining air in system to be purged out as described in the console operator manual.

Disconnecting Catheter from CoolGard 3000/Thermogard XP System

1. Stop circulation of saline through catheter.
2. Disconnect Start-Up Kit from catheter.
3. To maintain sterile connections, immediately cap off Luer connectors of both catheter and Start-Up Kit using sterile Luer caps or connect INFLOW and OUTFLOW Luers together.

Reconnecting Catheter to CoolGard 3000/Thermogard XP System

1. Remove Luer caps from Luer connectors of catheter and Start-Up Kit and discard or disconnect INFLOW and OUTFLOW Luers from each other.
2. Attach Start-Up Kit to Solex Catheter by connecting the male Luer of the Start-Up Kit to the female INFLOW Luer of the Solex Catheter and the female Luer of the Start-Up Kit to the male OUTFLOW Luer of the Solex Catheter. Assure that a sufficient amount of sterile saline is present at the ends of the hubs make an air free connection.
3. **WARNING: Failure to connect the Start-Up Kit correctly to the catheter could result in catheter failure.**
4. **WARNING: DO NOT use the INFLOW and OUTFLOW Luer fittings for standard central line infusion ports. They are for connection to the CoolGard 3000/Thermogard XP System ONLY.**
5. **Caution: Do not place any extra stopcocks in line that may be inadvertently shut off. This can cause line blockage and possible failure.**

Catheter Removal

1. Stop all pumping of saline through the catheter.
2. Disconnect Start-Up Kit from catheter. **Uncap or leave uncapped the saline INFLOW and OUTFLOW lumens of the catheter.** This will allow residual saline within the circuit to be expressed. As catheter is withdrawn, the balloons are compressed. **Saline within the balloons must be free to pass out of the balloon or the balloon will not deflate making the catheter difficult to remove.**
3. Place patient in supine position. Remove dressing. Remove sutures from suture site.
4. Attach a 20cc syringe to the catheter saline INFLOW hub. Pull and hold a vacuum for 15 seconds to allow residual saline to be removed from the catheter balloon section prior to starting to remove the catheter.
5. **WARNING: Do not place a luer cap on the OUTFLOW luer.**
6. Slowly remove catheter from patient. As catheter exits the site apply pressure with a dressing impermeable to air, e.g. Vaseline gauze.
7. **WARNING: Do not move catheter if resistance is felt. Check to ensure that the INFLOW and OUTFLOW lumens of the cooling circuit are NOT capped. If they are capped, uncap them and try removing the catheter again. If resistance is still encountered, an x-ray should be performed to identify the reason for the resistance.**

MRI INFORMATION

Non-clinical testing has demonstrated that the Solex Catheter is **MR Conditional**. Patients with the catheter inserted can be scanned safely under the following conditions

- Catheter must be disconnected from the CoolGard 3000 or Thermogard XP Console
- Static magnetic field strength of 1.5 T or 3 T only
- Maximum spatial gradient field of 720 Gauss/cm (7.2 T/m)
- Averaged whole body SAR of 2.9 W/kg for 15 minutes of scanning (per sequence)

MRI Artifacts

The presence of the Solex Catheter may cause moderate artifacts on the MRI image depending on the pulse sequence parameters used for MRI. However, the artifacts are confined to the position of the Solex Catheter and, as such, will not affect the diagnostic use of MRI imaging for head-only MRI examinations.

WARNING: The ZOLL Coolgard 3000 and Thermogard XP Consoles are MR Unsafe. Do not use in the MR Suite.

ZOLL Circulation, Inc.
2000 Ringwood Ave.
San Jose, CA 95131 U.S.A.

+1 (408) 541-2140 (main)

+1 (408) 541-1030 (fax)

Solex® is covered by United States Patent numbers: 6,749,625; 6,796,995. Also covered by European Patent numbers: 1,589,915; 2,174,625 and Japanese Patent number 4,527,121.

© 2014

GuideWire Instructions for Use

Caution: Federal law (USA) restricts this device to sale and use by or on the order of a physician.

Note: This information applies only to the use of guidewires in the Seldinger technique of catheter placement in the vasculature.

Warnings

The supplied guidewire is designed for single use only. Do not resterilize or reuse. Do not reinsert, once removed from patient.

Should resistance occur during insertion or withdrawal, DO NOT continue to move the guidewire. Determine the cause under fluoroscopy and take action as needed.

Use extreme caution when moving a guidewire through a stent. Use of a guidewire in stented vessels creates additional patient risk.

Cautions

Avoid withdrawing guidewire through metal needles; guidewire may shear.


Because of the delicate and fragile nature of guidewires, extra care in handling must be taken. Avoid bending or kinking. Do not use damaged guidewires.

Sufficient guide-wire length must remain exposed to maintain firm grip on guidewire at all times.

Dispenser

Every guidewire is provided in a dispenser package. Remove the guidewire anti-migration clip before dispensing the guidewire. Remove the guidewire protective cap immediately prior to guidewire use. Prepare the guidewire prior to insertion. It is recommended that the dispenser be filled with heparinized solutions (e.g. saline or dextrose) to bathe the guidewire during insertion.

Preformed "J" guidewire will resume shape when removed from product dispenser.

1 = Guidewire protector cap


2 = Guidewire anti-migration clip

Inspection


Inspect the guidewire prior to use and discard should any deformities be present in the guidewire. Guidewire placement should be routinely monitored by x-ray or fluoroscopic procedure.

Technique


1. Puncture vessel.
2. Insert the guide wire into the needle hub and gently advance 5-10cm of the guidewire into the punctured vessel. Navigate the guidewire to the desired position.
3. **Caution: Avoid rough or overly vigorous manipulation of the guidewire to prevent damage to the guide or the vessel.**
4. Remove the needle from the guidewire.
5. Dilate tissue and vessel with the dilator using a slight rotary motion.
6. Remove dilator. (Vessel dilator intended for vascular dilation only.)
7. Introduce catheter by sliding it over the guidewire.
8. Remove guidewire.

Solex Intravascular Heat Exchange Catheter

Solex Intravascular Heat Exchange Catheter Manifold and Distal Configuration

Solex Intravascular Heat Exchange Catheter Serpentine Balloon and Cross Section of the Inflated Solex Serpentine Balloon.

*The serpentine balloon has a cross sectional area of 54mm², which is approximately equal to the cross sectional area of an 8mm OD standard cylindrical balloon.

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As superfícies de contato do Solex Catheter com o sangue são tratadas com revestimento de heparina Applause™. Applause é marca comercial da SurModics, Inc., registrada no Departamento de Marcas e Patentes dos EUA (U.S. Patent and Trademark Office).
Het oppervlak van de Solex-katheter dat in aanraking komt met bloed, is behandeld met Applause™ heparinecoating. Applause is een handelsmerk van SurModics, Inc., gedeponeerd bij het U.S. Patent and Trademark Office (Amerikaanse octrooi- en handelsmerkbureau).
Povrchy katétru Solex Catheter, které jsou ve styku s krví, jsou opatřeny heparinovým potahem Applause™ Heparin Coating. Applause je ochranná známka společnosti SurModics, Inc. registrovaná u amerického Úřadu pro patenty a ochranné známky (U.S. Patent and Trademark Office).
Površine Solex Catheter-a koje dolaze u kontakt s krvlju tretirane su heparinskim premazom Applause™. Applause je zaštićen znak tvrtke SurModics, Inc., registriran u Američkom uredu za zaštitu patenata i robnih znakova. (U.S. Patent and trademark office).