Surface Start-Up Kit Instructions for Use

Surface Start-Up Kit Instructions for Use (EN).....	1
Návod k použití Souprava Surface Start-Up Kit (CS).....	10
Surface Start-Up Kit Brugsanvisning (DA)	18
Surface Start-Up Kit Gebrauchsanweisung (DE).....	26
Οδηγίες χρήσης του Surface Start-Up Kit (EL).....	34
Modo de empleo del kit de inicio de almohadillas superficiales (ES)	42
Potilaspinnan aloitussarjan käyttöohjeet (FI)	50
Mode d'emploi du kit de démarrage des coussinets de surface (FR)	58
Petunjuk Penggunaan Kit Awal Permukaan (ID)	66
Istruzioni per l'uso del Surface Start-Up Kit (IT).....	74
Gebruiksaanwijzing voor de opstartkit voor de onderlegger (NL).....	82
Bruksanvisning for startsett for overflatebehandling (NO)	90
Instruções de utilização do kit de inicialização da superfície (PT).....	98
Bruksanvisning till Surface Start-Up Kit (SV)	106
Hướng dẫn Sử dụng Bộ Khởi động Bề mặt (VI).....	114

Device Description

The Surface Start-Up Kit is a tubing set that connects the Thermogard XP® console to the ZOLL surface pads. The Surface Start-Up Kit is a closed-loop disposable that is composed of a pre-attached water reservoir, extension tubing, a pump loop, a heat exchange coil and an air trap assembly. All of these items work with the console to supply cool water to the surface pad to regulate patient core temperature.

Indications for Use

Vest: Thermogard XP® Temperature Management System: Temperature reduction in adult patients where clinically indicated, e.g. in hyperthermic patients.

VitalTemp™ Universal Pad*: When used with Thermogard XP or Thermogard HQ console: Temperature reduction in patients 2 years and older where clinically indicated, e.g. in hyperthermic patients.

* Not available in all regions.

Intended Purpose

The Surface Start-Up Kit intends to lower patient temperature.

Intended User

Intended Users for the Surface Start-Up Kit include physicians, physician assistants, and nurses.

Clinical Benefits

The primary clinical benefit of the Surface Start-Up Kit, used as part of the Thermogard system, is to lower body temperature.

Contraindications

None known.

Warnings

- Read these instructions for use carefully before connecting the surface pad to a patient.
- Systemic hypothermia risks. Systemic hypothermia may cause cardiac arrhythmia, patient shivering, or other system or organ complications. Systemic hypothermia should only be utilized under the supervision of a qualified physician.
- Surface Start-Up Kits (non-sterile) are not compatible with intravascular catheter use (sterile).
- The patient's temperature and skin condition must be checked, especially areas in contact with the surface pad, every 20 minutes or as directed by a physician. Patients vary in degree of sensitivity to cold, heat and pressure. Patients at greatest risk are those unconscious or incapacitated, persons on prolonged therapy, diabetics, children, persons with sensitive skin areas or poor circulation. Focus attention on all bony prominences.
- When a Thermogard XP console is used with a Surface Start-Up Kit and surface pad, the Treatment Mode should set to Fever. The console does not warm the patient when Fever mode is selected. The Lo patient temperature alarm limit ensures that an alarm occurs should the patient stop regulating his/her own body temperature. Such patients cool to room temperature. This can occur when the patient dies or becomes comatose. Investigate all patient temperature alarms.

- The Surface Start-Up Kit is MR Unsafe and, thus, is not allowed in the MR system room. Therefore, the Surface Start-Up Kit must be disconnected from the surface pad prior to moving the patient into the MR system room.

Cautions

- U.S. federal law restricts the device to sale on or by the order of a physician.
- Surface Start-Up Kits are not sterile and are intended for single patient use. Do not reprocess or sterilize.
- Single-use device – do not reuse. ZOLL Surface Start-Up Kits are single-use devices and may not be reprocessed or reused. The cyclical stresses of the roller pump on these items cause fatigue failures. Do not use Surface Start-Up Kits beyond the labeled usage time. The products will fail.
- Use care when fitting the surface pad to the patient. When priming is complete, continue to flow water through the surface pad before fastening the surface pad to the patient. Once water is fully flowing through the surface pad, it may be secured to the patient.
- Failure to connect both connectors to the surface pad may result in overpressure and leaks in the Surface Start-Up Kit or surface pad.
- If you notice a depleted reservoir bag or an air trap alarm, do not replace the water prior to identifying the location of the leak and taking the appropriate mitigation. Check for leaks in the Surface Start-Up Kit and the surface pad. If a leak is found, replace as necessary. If the surface pad is replaced, the Surface Start-Up Kit reservoir bag may be refilled with additional water as required.

Operating Life

Caution. Surface Start-Up Kit lifetime is 10 days. The designed operating lifetime for these components is 10 days of continuous operation on a single patient. If a patient must be treated for a longer period, a new Surface Start-Up Kit must be installed in the console. Failure to adhere to this time limit may cause the device to fail.

Setup Instructions

1. Initiate Console

1.0 POWER ON THE CONSOLE

1.1 SELECT "PRE-COOL"

1.2 CHOOSE NEW PATIENT (Yes or No)

Choose No if continuing a prior treatment.

1.3 SELECT TARGET TEMPERATURE

1.4 SELECT "FEVER" MODE

Caution. Always use "Fever" mode when using Thermogard surface treatment.

2. Place Surface Pad

2.0 PLACE SURFACE PAD

Vest: Place the surface pad under the patient with the white side up. Do not secure the Velcro straps at this time.

VitalTemp Universal Pad*: Place the pad under, over, or wrapped around the patient. Put the white side against the patient's skin. Do not secure the Velcro straps at this time.

* Not available in all regions.

Note. The surface pad may be disconnected from the tubing while positioning it under the patient.

3. Set Up Tubing

3.0 FILL RESERVOIR BAG

- A. Add water to the reservoir bag and fasten the cap.
- B. Hang the reservoir bag on the console hook.

3.1 INSERT COIL

- A. Open the coolant well cap.
- B. Place the coil in the well.
- C. Close the cap, passing the tubes through the slit in cap.

3.2 PLACE AIR TRAP

- A. Temporarily place the air trap in the air trap chamber.

3.3 LOAD PUMP TUBE

- A. Open the roller pump lid.
- B. Fully insert the tubing disc into the slot on the right side of the pump raceway.
- C. Lift the handle on the pump rollers.
- D. Firmly rotate the pump handle counter-clockwise while feeding the tube into the raceway.
- E. Fully seat the tube into the raceway.
- F. Press the handle down onto the rollers until it presses into its detent.

3.4 CONNECT SURFACE PAD

- A. Fully connect both fluid connectors to the surface pad.

Caution. The tubing must be fully connected to the surface pad (inlet and outlet) before it can be primed.

3.5 PRIME TUBING

- A. Close the roller pump lid.
- B. Remove the air trap from the chamber and hold it upside down.
- C. Prime the air trap by pressing and holding the Prime switch.
- D. When the air trap is completely full of water, release the Prime button.

3.6 INSERT AIR TRAP

- A. Turn the air trap right side up and insert it into the air trap chamber.
- B. Close the console top cover while ensuring the tubing is not pinched.

3.7 CONNECT PROBE

Connect a core temperature probe (third party) to the front of the console.

4. Start Treatment

4.0 PRESS STANDBY/RUN

Press the Standby/Run button on the console to start treatment.

Caution. Ensure the Velcro straps are not secured.

4.1 SET LO ALARM

Press the Menu button on the console, then set the patient Lo alarm to 1°C below the target temperature.

Warning. In Fever mode, the console will not warm the patient. Setting the Lo patient temperature alarm ensures that an alarm occurs if the patient's temperature drops.

4.2 SECURE SURFACE PAD TO PATIENT

Wait for the surface pad to fill with water. This typically takes 3-5 minutes. A small amount of air in the surface pad is acceptable and will dissolve over time. Confirm the water in the reservoir bag is still above the minimum level.

Vest: Wrap the surface pad flaps over the abdomen or chest. Secure and adjust to fit the patient using the enclosed stretchable Velcro straps. The surface pad should be snug but not constraining.

VitalTemp Universal Pad: Place the pad under, over, or wrapped around the patient. Secure with Velcro straps as required.

4.3 CONFIRM FLUID FLOW

Confirm that fluid is flowing through the tubing and surface pad by observing the rotation of the inline flow indicator (pinwheel).

This may take several minutes. Check for kinks and tap on the flow indicator as needed.

4.4 DURING TREATMENT

Check the patient temperature and skin condition, especially areas in contact with the surface pad, every 20 minutes or as directed by a physician.

Disconnect from Console and Dispose of Used Components

Caution. Avoid contact with used components. Handle used components as medical waste.

To remove used components from the console and dispose of them properly, follow these steps in the indicated order.

1. Place the console in Standby
 - A. If after discarding the used disposables, they will be replaced and therapy will continue:
 - Leave the console in Standby.
 - Discard the used disposables as indicated in steps 2-14 below.
 - Setup and initiate new treatment according to the respective instructions.
 - B. If treatment has concluded:
 - Turn off the power switch.
 - Discard the used disposables as indicated in steps 2-14 below.
2. Disconnect the Surface Start-Up Kit connectors from the surface pad.
3. Position a large empty medical waste container or collection bag near the console.

Caution. Do not disassemble Surface Start-Up Kit. Do not disconnect the tubing connecting the components. They are removed and disposed of as a complete unit.
4. Open the top cover of the console. Open the transparent top cover of the roller pump.
5. Remove the lid from the coolant well and set it aside in a clean location.
6. Lift the handle on the pump rollers.
7. Grasp the pump tubing and gently pull it up and out of the channel while rotating the pump head.
8. Pull the tubing out of the pump head. Press the handle down onto the rollers until it press into the detent and close the top cover of the roller pump.
9. Lift the heat exchange coil out of the coolant well. Hold the coil above the coolant well to allow coolant to drip back into the well.
10. Pull the air trap up and out of the chamber.
11. Remove the water reservoir from the hook.
12. Bundle up the surface pad and Surface Start-Up Kit components, still connected together, and gently deposit them in the medical waste container.
13. Use a tissue to wipe up any spilled coolant from the top of the coolant well. Place the tissue in the medical waste container.
14. Replace the insulated lid on the coolant well.

Removal and disposal of the used components is complete. The console may now be stored or moved to its next treatment location.

Troubleshooting

- If pad surface is not cold after approximately 5 minutes, check for kinks and other restrictions to water flow in the tubing and the surface pad.
- It is normal for some amount of air to release from the surface pad during the procedure. This air will be captured in the water reservoir.
- Refer to the surface pad Instructions for Use for full troubleshooting guidance.

Compatibility

- The Surface Start-Up Kit works in conjunction with the Thermogard XP console.
- The Surface Start-Up Kit connects to the surface pad.

Storage

- Store the Surface Start-Up Kit in a cool, dry environment.

Serious Incident Reporting

- If a serious incident related to the device occurs, immediately report the incident to a healthcare professional.
- For healthcare professionals, immediately report any serious incident to the applicable competent authority and ZOLL.

Symbols

	Warning/Caution
	Do Not Re-use; Single Use
	Non-Sterile
	Lot or Batch Number
	Catalog Number
RX ONLY	U.S. Federal law restricts this device to sale by or on the order of a physician
	Conformité Européenne
	Unique Device Identifier
	Importer
	Authorized representative for Switzerland
	MR Unsafe

Description	Catalog Number
Thermogard XP Surface Start-Up Kit	8700-000936-01

Copyright ©2025 ZOLL Circulation. All rights reserved. ZOLL, Thermogard, Thermogard XP, and VitalTemp are trademarks or registered trademarks of ZOLL Circulation and/or its subsidiaries in the United States and/or other countries. All trademarks are the property of their respective owners.

Patent: www.zoll.com/patents

	Consult Instruction for Use
	Manufacturer
	Date of Manufacture
	Do not use if package is damaged and consult instructions for use
	Keep Away from Heat/Sunlight
	Keep Dry
	Quantity
	Medical Device
	Use by date
	Authorized representative in the European Community/European Union

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Copyright ©2025 ZOLL Circulation. Všechna práva vyhrazena. ZOLL, Thermogard, Thermogard XP a VitalTemp jsou ochranné známky nebo registrované ochranné známky společnosti ZOLL Circulation a/nebo jejích dceřiných společností ve Spojených státech amerických a/nebo v jiných zemích. Všechny ochranné známky jsou majetkem příslušných vlastníků.
Copyright ©2025 ZOLL Circulation. Alle rettigheder forbeholdes. ZOLL, Thermogard, Thermogard XP og VitalTemp er varemærker eller registrerede varemærker tilhørende ZOLL Circulation og/eller virksomhedens affilierede selskaber i USA og/eller andre lande. Alle varemærker tilhører deres respektive ejere.
Copyright © 2025 ZOLL Circulation. Alle Rechte vorbehalten. ZOLL, Thermogard, Thermogard XP und VitalTemp sind Marken oder eingetragene Marken von ZOLL Circulation und/oder ihren Tochtergesellschaften in den USA und/oder anderen Ländern. Alle Marken sind Eigentum ihrer jeweiligen Inhaber.
Copyright ©2025 ZOLL Circulation. Reservados todos los derechos. ZOLL, Thermogard, Thermogard XP y VitalTemp son marcas comerciales o marcas registradas de ZOLL Circulation o sus entidades subsidiarias en EE. UU. o en otros países. Todas las marcas comerciales son propiedad de sus respectivos titulares.
Copyright © 2025 ZOLL Circulation. Tous droits réservés. ZOLL, Thermogard, Thermogard XP et VitalTemp sont des marques commerciales ou déposées de ZOLL Circulation et/ou de ses filiales aux États-Unis et/ou dans d'autres pays. Toutes les marques sont la propriété de leurs titulaires respectifs.
Copyright ©2025 ZOLL Circulation. Tutti i diritti riservati.