

# When your response team needs more than an AED, it's time to consider the R Series



The ZOLL® R Series® monitor/defibrillator provides clinicians with comprehensive support for resuscitation. This includes cutting-edge technology to help you meet current Guidelines for achieving high-quality CPR, as well as OneStep™ electrodes that simplify and speed therapy. And because the R Series conducts a comprehensive self-test daily, you can be assured that it will always be Code-Ready®.

## **CPR Feedback—The New Standard**

Every victim of cardiac arrest requires CPR—high-quality CPR. The R Series has Real CPR Help® technology to provide coaching on CPR depth and rate, so you know you're meeting the 2015 American Heart Association (AHA) Guidelines. With CPR Dashboard™, detailed visual feedback, including the depth and rate of each compression is displayed on the R Series screen.

## **One-Step Electrodes—Simple and Comprehensive**

Monitor, pace, defibrillate, and get Real CPR Help for all patients—pediatric and adult—using only one type of electrode. Electrodes are preconnected and are automatically tested for an expiration date.

## **Comprehensive Readiness Checks**

Problems during a code are often related to compromised supplies and batteries. The R Series extends testing far beyond a basic test shock to check more than 110 individual indicators of readiness, including batteries, cables, and electrodes.

**ZOLL**



### R Series ALS

Basic Defibrillation with OneStep Pacing  
ZOLL Part Number – 35020000001110012

The 2015 AHA Guidelines emphasize the importance of minimizing interruptions in CPR. In addition to Real CPR Help, the R Series ALS provides ZOLL's unique See-Thru CPR® technology. See-Thru CPR reduces the duration of interruptions by allowing clinicians to see whether an organized shockable rhythm has developed during compressions, which reduces the duration of pauses in CPR.



### R Series BLS

Basic Defibrillation with OneStep Pacing  
ZOLL Part Number – 35020000001130013

When the R Series BLS powers on, it shows first responders a simple, one-button AED interface. If an ALS-trained rescuer arrives, with the press of a button, the BLS unit turns into a full-function manual defibrillator. All user-chosen advanced monitoring options are functional and, if installed, the pacing controls are revealed. This entry-level model allows medical facilities to deploy R Series technology.

#### The R Series ALS and R Series BLS includes:

- 1 SurePower™ rechargeable Li-ion battery pack
- 2 sets of OneStep Complete electrodes
- Real CPR Help
- Code readiness testing system
- ZOLL Rectilinear Biphasic™ waveform technology
- Compact flash card slot
- Integrated AC power (120 VAC/60 HZ)
- 5-year warranty
- On-site in-service included

#### Additional parameters available:

- Pulse oximetry (SpO<sub>2</sub>)
- Non-invasive blood pressure (NIBP)
- Capnography (EtCO<sub>2</sub>)
- Expansion pack (See-Thru CPR, CPR Dashboard, Wi-Fi capability)

To find out more about how ZOLL can help, contact your ZOLL representative. If you aren't sure who that is, call **888-851-8584** or email [distributorsupport@zoll.com](mailto:distributorsupport@zoll.com) to find out.