

Real CPR Help Travel Trainer

Experience Real CPR Help with Real-Time Feedback



The Real CPR Help® Travel Trainer combines a simulated torso, **CPR-D•padz**®, and simulator in an easy-to-use, compact package ideal for demonstrating or practicing CPR with Real CPR Help.

ZOLL
Advancing Resuscitation. Today.®

Real CPR Help Travel Trainer includes:

- One-piece zippered carry-case
- Simulated **CPR-D•padz**[®] and embedded heart rhythm simulator
- 9-volt battery
- 6-minute in-service video (English only)

ZOLL Medical Corporation Worldwide Headquarters

269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-348-9011

ZOLL International Offices

Australia
+61 2 8424 8700
www.zoll.com.au

Austria
+43 650 413 6222
www.zollmedical.at

Canada
905-629-5005
www.zoll.com

China
+86 21 6226 9505
www.zoll.com

France
+33 (1) 30 05 14 98
www.zoll.fr

Germany
+49 2236 87 87 0
www.zollmedical.de

Middle East and Africa,
+96 278 800 4411
www.zoll.com

The Netherlands
EU Authorized Representative
+31 (0) 488 41 11 83
www.zoll.nl

Russia
+7495 936 2338
www.zoll.com

United Kingdom
+44 (0) 1928 595 160
www.zoll.com

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.

No Manikin or Simulator Needed

Everything needed to operate the Real CPR Help[®] Travel Trainer is included in this one-piece accessory. The simulator—which is securely embedded in the thick vinyl of the Travel Trainer—is manufactured by Symbio, the leading provider of cardiac heart rhythm simulators. Simply connect the Real CPR Help Travel Trainer to a ZOLL[®] clinical AED Plus[®] or AED Pro[®] (sold separately) to experience ZOLL's Real CPR Help with real-time feedback for each chest compression.



One-piece zippered carry-case simplifies transport and storage.

Real CPR Help Travel Trainer Specifications

Dimensions

In Travel Trainer Carry Case: 9 ½" H x 18" W x 3" D (24.1cm H x 45.7cm W x 7.6cm D)

Flat: 32" H x 17¾" W x 3" D (81.3cm H x 45cm W x 7.6cm D)

Weight: 1 lb 12 oz (0.79 kg)

Battery: 9V battery with standard PP3 connectors (included)

Simulated Heart Rhythms

- No impedance - electrode pads not attached
- Asystole - electrode pads attached, no shock advised
- Ventricular Fibrillation - electrode pads attached, shock advised

Compatibility

- ZOLL clinical AED Plus
- ZOLL clinical AED Pro



Specifications subject to change without notice.

©2007 ZOLL Medical Corporation. All rights reserved. "Advancing Resuscitation. Today.™", AED Plus, AED Pro, **CPR-D•padz**, Real CPR Help, and ZOLL are registered trademarks of ZOLL Medical Corporation. Printed in U.S.A. 120710 9656-0203

ZOLL[®]
Advancing Resuscitation. Today.™