

# RescueNet® DeviceDashboard

## Maintenance and support simplified

RescueNet DeviceDashboard asset management software helps you stay on top of defibrillator maintenance across multiple devices and campuses. This cloud-based software saves time, adds value, and helps ensure defibrillators are ready the second they're needed.



RescueNet DeviceDashboard asset management software is hosted securely by ZOLL® Online, giving you immediate access to device readiness reports using your unique credentials.

# ZOLL®

# RescueNet DeviceDashboard: Ensuring Code Readiness



## Robust readiness checks:

- Defibrillation, pacing, and ECG modules
- AC power connection
- Battery status
- 30-joule shock test
- Electrode attached and within expiration date
- Automatic daily clock synchronization



## Data securely hosted by ZOLL® Online:

- No need for a separate server, on-premises maintenance, or hardware updates
- Product updates, security enhancements, and bug fixes installed automatically
- Streamlined access to ZOLL technical support



## Hospital-wide views:

- Comprehensive daily report on fleet's state of readiness
- Highly customizable
- History of readiness logs to comply with The Joint Commission
- Secure login from any computer or mobile device with internet access

## Daily Status Update

The R Series® and X Series® monitor/defibrillator and the ZOLL AED 3® BLS defibrillator conduct comprehensive self-tests of various individual indicators of readiness. Using standard 802.11 a/b/g/n wireless technology, they securely transmit this information to DeviceDashboard. This data is compiled into at-a-glance reports that can be viewed using Microsoft Edge, Google Chrome, or Apple Safari on a PC or laptop.

The screenshot shows the ZOLL DeviceDashboard interface on a laptop. The interface includes a sidebar with filters for Device Type, Site ID, and Assigned Location. The main content area displays a table of devices with columns for Device, Assigned Location, Last Check-in, Ready for Use, and Notifications. Callouts point to various features:

- Location of entire fleet across a hospital
- Last defibrillator self-test
- Self-test results: Ready, Not Ready, or Fail
- Electrode expiration date
- Exportable .csv asset reports
- Type of device:
  - R Series®
  - X Series®
  - ZOLL AED 3® BLS
- Consolidated, customizable email alerts and status notifications

Find out more at [zoll.com/products/data/hospital](https://zoll.com/products/data/hospital)