
Check your AEDs from anywhere

With the CheckAED® app, you can quickly and efficiently perform and input routine inspections from anywhere. Simply download the app to your smartphone or tablet. When it's time to check each AED, inspect its status indicator and log it in PlusTrac with a single touch. If anything needs correction or replacement, your program manager will be notified right away.





Rescue Ready® Services Better Management for Better Outcomes

To keep your AED program running smoothly, you need the right partner and services. Keeping your volunteers ready, devices maintained, AED program compliant and, above all, your people safe each and every day can be easy with PlusTrac and Rescue Ready Services.

PlusTrac leading-edge software and flexible service offerings are here to support you.



To learn more, visit us at www.zoll.com/plustrac.

ZOLL MEDICAL CORPORATION

269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | www.zoll.com

Copyright © 2024 ZOLL Medical Corporation. All rights reserved. CheckAED, ICPR, PlusTrac, Powerheart, Rescue Ready, ZOLL AED 3, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

¹Drezner JA, et al. Br J Sports Med. 2013;47:1179-1183

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.











Increasing a sudden cardiac arrest victim's chance of survival requires more than just placing an AED on the wall.

- Devices must be **ready.**
- Responders must be **ready**.
- Programs must be **compliant.**

PlusTrac Keeps You on Track

PlusTrac is an interactive, web-based solution that enables seamless, reliable AED management for programs of all sizes. From registering your AED and documenting your responder training to recall notifications, legislative updates, and AED prescriptions, PlusTrac helps streamline administration and provides visibility to all facets of your program.

Whatever brand of AED you've chosen, PlusTrac will help you manage your whole program.



Studies show that well-implemented AED programs can significantly improve survival rates.¹

Know You're Ready - At a Glance

PlusTrac makes it easy to manage and track your AEDs, their pad and battery status, inspections, and responder training and certifications. And customized reporting facilitates teamwide communication.

Roles for Collaborative Organizations

The bigger the organization, the greater the need to set clear roles and responsibilities for team members. PlusTrac provides views that streamline information according to individual roles, creating a platform that scales for even the largest teams.

Software and Service for the Win


Effective program management requires more than software. From FDA-required AED prescriptions and medical oversight to EMS agency notification and post-event support, PlusTrac delivers everything you need for worry-free AED program compliance.

Connected Devices Build Confidence

WiFi-connected ZOLL AED 3 devices self-report AED status, battery use, and electrode expiration, along with post-event data that PlusTrac translates into usable, actionable alerts. Similarly, ZOLL Powerheart G5 devices connected via the Rescue Ready Monitor transmit similar data.








A Scalable Solution for Organizations of All Sizes – from 1 AED to 1,000 Regardless of how many AEDs or rescuers are in

Regardless of how many AEDs or rescuers are in your program, managing them does not need to be complex. PlusTrac can help you with a single AED in a single facility or hundreds of AEDs in multiple facilities with thousands of rescuers to train.





MEDICAL DIRECTION

- State and local compliance guidance
- AED maintenance intervals
- Post-use medical reviews



AUTOMATIC ALERTS

- Pad and battery expiry notices
- Responder certification expirations
- Manufacturer recalls (any device)



REAL-TIME INFO

Remote WiFi reporting for connected ZOLL AED 3® devices and Powerheart® G5 AEDs with Rescue Ready® Monitor



USER-FRIENDLY DASHBOARD

- Interactive management console
- Easy access to program status
- High-level overview with drilldown to more detailed information