

ECG Electrodes

Ag/AgCl Pregelated Monitoring Electrodes

Contains
3

AW-R2246-03 Rev. B



ZOLL

ZOLL Medical Corporation
269 Mill Road
Chelmsford, MA 01824-4105 USA
900-948-9011 Made in USA

EC|REP ZOLL International Holding B.V.
Newtonweg 18 6662 PV ELST
The Netherlands Tel. +31 (0) 481 366410



- 2 Cautions:
- Do not use if gel is dry.
 - Replacement of electrodes after each 24 hour period is suggested. Vary placement to minimize skin sensitivity.

