

# ECG Electrodes

Ag/AgCl Pregelled Monitoring Electrodes

Contains  
**5**

AW-R2246-05 Rev. A



1



2



3



4

## ZOLL

ZOLL Medical Corporation  
269 Mill Road  
Chelmsford, MA 01824-4105 USA  
900-348-9011 Made in USA

EC REP ZOLL International Holding B.V.  
Newtonweg 18 6662 PV ELST  
The Netherlands Tel. +31 (0) 481 366410

### Cautions:

- Do not use if gel is dry.
- Replacement of electrodes after each 24 hour period is suggested. Vary placement to minimize skin sensitivity.

