


# Stat-padz® II

## Training Electrodes

### Instructions for Use

1. Attach **Stat-padz® II** Training Electrodes to AED Plus® Trainer or ZOLL® training cable.
2. Attach the **Stat-padz® II** Training Electrodes to any CPR training manikin, pressing the electrodes firmly to the surface and moving any air pockets to the outer edges.



### IMPORTANT:

- Not for clinical use.
- For training purposes only.
- Use with AED Plus® Trainer or ZOLL® training cable.

OPERATING TEMPERATURE: 0°C to 50°C (32°F to 122°F)  
SHORT TERM STORAGE TEMPERATURE: -30°C to 65°C (-22°F to 149°F)  
LONG TERM STORAGE TEMPERATURE: 0°C to 35°C (32°F to 95°F)



NOT MADE WITH  
NATURAL RUBBER LATEX.

R2001-38 Rev. J

# ZOLL®



ZOLL Medical Corporation  
269 Mill Road  
Chelmsford, MA 01824-4105 USA  
800-348-9011 Made in USA



ZOLL International Holding B.V.  
Newtonweg 1B 6662 PV ELST  
The Netherlands Tel. +31 (0) 481 366410