

Product Bulletin

4 December 2025


RescueNet® CodeWriter - Rapid Response Release


Overview

We are pleased to announce the release of **CodeWriter version 1.0.15**. This update introduces significant enhancements to improve usability and streamline workflows. Key highlights include:

- **Rapid Response Documentation:** Optimized to support both resuscitation and rapid response events.
- **User Interface Enhancements:** Improved layout and navigation for a more intuitive experience.
- **Case Upload Optimization:** Faster and more reliable submission process.

New Features


- **Event Details:** Users can now enter supplemental event information before completing active documentation. To access these fields, tap the new  icon located in the top-right corner of the home screen.



Product Bulletin

4 December 2025


Expanded Vital Signs Documentation: Clinicians can now record additional fields within the vital sign screen including - cardiac rhythm, O₂ source, flow rate, and FiO₂.



- **Expanded Cardiac Rhythm Selection:**

The rhythm list now includes:

- 1st Degree Block
- 2nd Degree Block (Type 1)
- 2nd Degree Block (Type 2)
- Atrial Fibrillation with RVR




Product Bulletin

4 December 2025


- **New Procedure Options**

- Reorganized and added additional airway procedures and O₂ delivery sources within the procedure submenus.
- Added CPR Backboard and ECMO procedures to the 'other procedures' list.



- **New Medication Options**

- Added: Adenosine, Dopamine, Etomidate, Ketamine, Norepinephrine, Propofol, and Succinylcholine.
- New IV Fluids Rapid Infusion Choices: Pressure Bag and Rapid Volume Infuser.




Product Bulletin

4 December 2025

- New 'Complete and Review' Workflow**

A redesigned workflow is now available for completing and reviewing documentation.


Important: Once a user completes this workflow, they cannot return to the primary documentation screen.



* The appearance of the user login screen (3rd image) may vary based on your institution's user management system.

- Automatic Upload Retry:**


- When a user taps 'Upload Event' on the 'Checklist' screen, the app will now automatically retry the upload if a network issue occurs.
 - Retries: 4 attempts with 15 seconds between attempts
 - If all attempts fail, the file moves to the Pending Uploads list.




Product Bulletin

4 December 2025

- **Manual Upload Retry:** Users can now manually upload files from the 'Pending Uploads' folder. Tapping 'Upload Events' will attempt to upload all pending event files.



- **Review Event Log Workflow Enhancements:** Users can now add items directly from the 'review event log' screen while reviewing documentation. This workflow mirrors the process on the primary documentation page.



Product Bulletin

4 December 2025




- **Full Feature List:** For a complete list of enhancements in version 1.0.15, please refer to the release notes available in the iOS and Android app stores.

Technical Details

- **Availability**
CodeWriter 1.0.15 is now available for iOS and Android devices in their respective app stores.
- **Compatibility**
This version is compatible with On-premise CodeNet version 10.2 or newer.
- *****End of Support Notice*****
On-premise CodeNet will reach **end of support on December 31, 2025**. If you have not yet migrated to the cloud, please contact your local territory manager or ZOLL Support to schedule your migration.

Contact

ZOLL Support teams are available to assist with any questions and to ensure a smooth implementation of these changes. If you need assistance, please contact our support team.

ZOLL Customers Within the US		1-800-348-9011		supportdata@zoll.com
ZOLL Customer Outside the US	Contact your local ZOLL Service Center or Authorized Distribution Partner			
		If you need assistance, please email ZOLL Service		intlservice@zoll.com