

# Automated External Defibrillators

# ZOLL®




When they need CPR,  
you need ZOLL.

# Confidence that comes only from **Real CPR Help.**

## **AED Plus®**

Step-by-step support for  
the entire Chain of Survival

*The ZOLL AED Plus gives even infrequent rescuers the confidence they need to help save lives. Simple audio prompts and illustrations reinforce every step in the resuscitation process.*




*The one-piece adult electrodes (CPR-D-padz®) with Real CPR Help® technology make pad placement fast and accurate. Real CPR Help can more than double a victim's chances of survival.<sup>1</sup>*

## ***"Push Harder." "Good Compressions." "Press Shock Button."***

Choosing ZOLL® for your next automated external defibrillator (AED) gives your patients yet another reason to trust you. As the brand hospitals choose most for professional defibrillators, ZOLL has earned the respect of users at every skill level. That's because our AEDs do more than provide the most advanced defibrillation technology.

The AED Plus® and the AED Pro® are the only automated external defibrillators to provide Real CPR Help®, real-time feedback on the depth and rate of chest compressions during CPR. In addition to audio prompts, both the AED Plus and AED Pro feature on-screen instructions, and a bar gauge that shows how deep to push in order to optimize compression depth. Helping rescuers to provide high-quality CPR is key because it has become increasingly clear just how important CPR is to survival from sudden cardiac arrest.



### **AED Pro®**

*A rugged, advanced AED  
for professional users*



Chest compression  
depth gauge



*The ZOLL AED Pro, whether used in a medical facility or in the field, is designed to perform in any environment, and is used by all types of rescuers. With Real CPR Help and See-Thru CPR®, it provides the advanced functionality that professional rescuers and services require.*





### AED Plus

- Real CPR Help provides real-time feedback on the depth and rate of CPR chest compressions
- Supports the entire Chain of Survival, with easy-to-understand audio prompts and illustrations
- Has lowest total cost of ownership, once installed, because CPR-D-padz and batteries last five years, if not used
- Rugged design that resists dust and direct water spray (IP55 rating)
- Powered by 10 lithium batteries available from retail stores everywhere



### AED Pro

- Real CPR Help provides real-time feedback on the depth and rate of CPR chest compressions
- See-Thru CPR filters CPR artifact to display patient's underlying ECG rhythm
- Advanced capabilities, including 3-lead patient monitoring and manual override
- The only AED that has an IP55 rating for dust and water resistance and passes the 1.5-meter drop test
- Battery and electrode compatibility with ZOLL's line of professional defibrillators

<sup>1</sup>Bobrow BJ, et al. *Ann Emerg Med*. 2013 Jul;62(1):47–56.e1. Epub 2013 Mar 7.

## ZOLL MEDICAL CORPORATION

An Asahi Kasei Group Company | 269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | [www.zoll.com](http://www.zoll.com)

Copyright © 2016 ZOLL Medical Corporation. All rights reserved. AED Plus, AED Pro, Real CPR Help, See-Thru CPR, CPR-D-padz, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

Printed in U.S.A.  
MCN PP 1606 0184

For subsidiary addresses and fax numbers, as well as other global locations, please go to [www.zoll.com/contacts](http://www.zoll.com/contacts).

**ZOLL**®