


Electrode Placement Guide—ZOLL M Series


Cardiac Emergencies

POSITION: ANTERIOR (APEX-LATERAL) / ANTERIOR (STERNUM)



- Emergency Defibrillation
- Emergency Ventricular Dysrhythmia Cardioversion


POSITION: ANTERIOR (APEX-FRONT) / POSTERIOR



- Emergency Defibrillation
- Emergency Ventricular Dysrhythmia Cardioversion
- Pacing

Elective Cardioversion

POSITION: ANTERIOR (STERNUM) / POSTERIOR



- Atrial Dysrhythmia Cardioversion