

solid gel **pedi-padz**™

for Emergency Situations



ZOLL
It's about time.

Expressly Designed for ZOLL Systems.

Our electrodes are specifically designed, tested, and validated for use with ZOLL defibrillators and pacemakers.

ZOLL systems offer complete capability. With ZOLL electrodes and standard cable connectors, you are assured of a complete system. No need for adapters, special connectors, or modified cables. There are no compromises with ZOLL electrodes.

Electrodes are also a key component in ZOLL's superior pacing technology, achieving higher capture rates at lower mean thresholds.

ZOLL—the name you can rely on for high performance, proven efficacy, and ease-of-use.

ZOLL Medical Corporation

269 Mill Road
Chelmsford, Massachusetts 01824-4105 U.S.A.
800-348-9011
978-421-9655
978-421-0025 Telefax
www.zoll.com

ZOLL Medical European Operations

Dodewaard, The Netherlands
+31 488 4111 83
+31 488 4111 87 Telefax

ZOLL Medical U.K. Ltd.

Cheshire, England
+44 1925 846 400
+44 1925 846 401 Telefax

ZOLL Medical Germany

Cologne, Germany
+49 221 3989 340
+49 221 3989 336 Telefax

ZOLL Medical Canada

Mississauga, Ontario, Canada
905-629-5005
905-629-0575 Telefax

ZOLL Medical France

Coignières, France
+33 1 30 05 1495
+33 1 30 49 0693 Telefax

ZOLL Medical Australia Pty. Ltd. and Far East Regional Office

New South Wales, Australia
+61 2 9420 8733
+61 2 9420 9834 Telefax

ZOLL Medical Middle East and Africa Region Office

Athens, Greece
+30 210 62 36691
+30 210 62 36692 Telefax

ZOLL Medical Japan

Tokyo, Japan
+81 45 905 2864
+81 45 905 2864 Telefax

ZOLL Medical Latin America

Parkland, FL 33067
954-345-4224
954-345-2684 Telefax

www.zoll.com

Meeting the challenges of pediatric patients and the rugged demands of the EMS environment.

pedi-padz Multi-Function Electrodes are the most technologically advanced Pediatric Multi-Function Electrodes designed to meet all Pediatric ACLS criteria. The high viscosity "solid" gel provides excellent skin coupling, lowers resistance quickly, virtually eliminates dermal damage, and has no storage limitations (does not have to be stored flat) making it an excellent choice for EMS applications.

pedi-padz Multi-Function Electrodes offer anatomically correct diagrams on the packaging to facilitate proper pad placement.

Solid gel **pedi-padz** Specifications

Shelf Life:	24 months
Conductive Gel:	Polymer hydrogel
Conductive Element:	Tin
Packaging:	Paper/poly/foil laminate
Impedance Class:	Low
Length, Lead Wires:	18.0 in (45.7 cm)
Pacing Time:	1 hour

FRONT (Apex) Electrode

Shape:	Round
Diameter, Total:	3.8 in (9.7 cm)
Area, Total:	10.8 sq in (69.3 sq cm)
Diameter, Conductive Gel:	2.9 in (7.2 cm)
Area, Conductive Gel:	6.4 sq in (41.1 sq cm)

BACK (Sternum) Electrode

Shape:	Rectangular
Length, Total:	5.4 in (13.8 cm)
Width, Total:	3.5 in (8.9 cm)
Area, Total:	19.0 sq in (122.6 sq cm)
Length, Conductive Gel:	3.6 in (9.1 cm)
Width, Conductive Gel:	2.2 in (5.5 cm)
Area, Conductive Gel:	7.8 sq in (50.3 sq cm)

Case Part Number:	8900-3000-01
Case Count:	6 pairs

pedi-padz have a backing material of closed cell polyethylene foam and connect via a pre-attached lead with keyed ZOLL connector.

Related Product:	Pediatric ECGs
Case Part Number:	8900-1003-01
Electrodes/pouch:	3
Pouches/case:	100



pedi-padz—specially sized to meet AHA pediatric ACLS criteria.

ZOLL

It's about time.