

CheckAED Mobile App

Need an easy and convenient way to manage your AED program? PlusTrac® program management includes the CheckAED® mobile app that lets you quickly and efficiently log routine inspections from any location.



With the CheckAED mobile app, you can easily:

- Inspect each of your AEDs on the spot
- Scan AED QR codes to speed inspection logging
- Receive alerts and notifications on your devices
- View your personalized dashboard
- Access customer support

Simplify inspections with QR code scanning

It's easy. From your web-based PlusTrac program management dashboard, print QR code labels for each of the AEDs in your fleet. Simply attach each QR code label to the appropriate AED, cabinet, or carry case.

When it's time to inspect your AEDs, use the CheckAED app on your mobile device to scan each device's QR code. The app matches the scanned serial number to an AED in your account, allowing you to easily log inspection details.

Ongoing management and inspections ensure your AEDs are ready and compliant.

Current customers can begin using the app right away. Simply download and enter your PlusTrac login credentials. Not yet a customer? Register for a **free trial**.

Give the CheckAED app a try to see how easy AED inspections can be.

*Contact your local representative for account access

ZOLL Medical Corporation | Chelmsford, MA, USA | 800-348-9011 | www.zoll.com

Copyright © 2022 ZOLL Medical Corporation. All rights reserved. CheckAED, PlusTrac, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries.
MCN PP 2205 0681

ZOLL®

